

CLEAN ROOM WASTE CONTAINER

CLEAN ROOM WASTE CONTAINER



Stainless steel waste container for disposing of used radioactive sources and/or needles.

The top is removable for the replacement of the disposable orpagard container. The top is fitted with a flip open lid to provide easy access to the opening in the container.

The lid can be operated with the back of the hand to avoid possible contamination. The waste container is shielded with lead which is encapsulated in a seamless stainless steel body, this also allows effective sterilisation idea for a clean room environment.

The stainless steel waste container comes in either 6mm Pb or 10mm Pb.

Code	Description	Lead Thickness (mm)	Internal dimensions LxWxH (mm)	External dimensions LxWxH (mm)	Approximate weight (kg)
6BT-CRWC6-OP2	For Orpagard 2 Litre	6	170x170x229	245x220x320	17
6BT-CRWC10-OP2	For Orpagard 2 Litre	10	170x170x229	245x220x320	21

CUSTOM SIZES ARE AVAILABLE ON REQUEST