

TUNGSTEN SYRINGE SHIELDS

UNI-LOCK® I TUNGSTEN SYRINGE SHIELDS

Original Uni-Lock® syringe shield.

Code	Weight (g)	Length (mm)	Int Ø (mm)	Ext Ø (mm)
1BT-U101A	120	72	8.8	12.4
1BT-U101B	140	77	8.8	12.4
1BT-U101L	165	75	11.5	16
1BT-U102A	140	53	11.5	16
1BT-U102B	145	59	11.5	16
1BT-U103A	150	65	11.5	16
1BT-U103B	160	68	11.5	16
1BT-U105A	190	60	15.4	20
1BT-U105B	195	65	15.4	20
1BT-U105C	200	68	15.4	20
1BT-U110A	250	75	18.2	22.5
1BT-U110B	260	85	18.2	22.5
1BT-U120	285	95	23	27

CUSTOM SIZES ARE AVAILABLE ON REQUEST

Note: If the make of syringe is not specified on an order, the longest length of that particular size will be supplied for radiation safety reasons.

