

TRANSPORT BAG



PRO-TEC® PIG TRANSPORT BAG

The durable, nylon, waterproof shipping bag is designed to accommodate up to 11 Pro-Tec® Pigs containing syringes. It is the connecting piece in the Biodex solution, ensuring safe, convenient handling of syringes from the pharmacy and back again. The system meets DOT 7A Type A packaging requirements.

Dimensions: 10.5x5.75x11.5" HxWxD (267x16x292mm)

Material:

- Exterior: Nylon
- Interior: Polyester



Code	Description
1001-754	Shipping Bag, Pro-Tec® Pig